



## SAFE LIFT: STOP FOR SAFETY



**Before placing resident in sling**



Get **SET**

**S**ling type and size is correct

**E**quipment is ready to receive resident

**T**est Battery

---

*After resident is in sling and loops attached to lift:*

*Lift approx. 2 inches and **S.T.O.P.***



**STOP**

**S**ling type and size is correct (*one last check*)

**T**he loops are fully secured

**O**pened base with wheels locked (if applicable)

**P**ositioning is correct and comfortable



**LIFT**

Staff & Resident (if able) give the green light to **LIFT**