

SAFE LIFT: STOP FOR SAFETY ORGANIZATION PLEDGE



By taking this pledge, I commit to implement the Safe Lift: STOP for Safety process in my organization to support our staff and keep our residents safe from harm.

Sign the
pledge by
using the QR
code.



Before placing resident in sling



Get SET

- S**ling type and size is correct
- E**quipment is ready to receive resident
- T**est Battery

*After resident is in sling and loops attached to lift:
Lift approx. 2 inches and **S.T.O.P.***



STOP

- S**ling type and size is correct (one last check)
- T**he loops are fully secured
- O**pened base with wheels locked (if applicable)
- P**ositioning is correct and comfortable



LIFT

Staff & Resident (if able) give the green light to **LIFT**